



A11.1020 LED 400x Achromatic Student Binocular Biological Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1020
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Eyepiece: WF10x/16mm
- Nosepiece: Triple Nosepiece
- Objective: Achromatic Objective
- Stage: Plain Stage With Slide Clip, Size 95*95mm
- Focusing: Coaxial Coarse & Fine Adjustment
- Condenser: N.A. 0.65 With Disc Diaphragm
- Illumination: Upper Flexible LED, Bottom LED Brightness Adjustable
- Power: Power By AC 110-240V Or DC 3xAAA Batteris (Batteries Not Included)
- Highlight: **400x Achromatic Student Binocular Biological Microscope**
, LED Student Binocular Biological Microscope,
Achromatic Student Binocular Biological Microscope

OPTO-EDU



More Images

OPTO-EDU



OPTO-EDU



Product Description

Compact Design Full Metal Frame Student Biological Microscope
Monocular, Binocular Head Available, With WF10x/16mm, WF25x/9mm Eyepiece
3 Holes Nosepiece With Achromatic Objective 4x,10x,40x(S)
Plain Stage 95*95mm, Coaxial Coarse & Fine Focusing
Up & Bottom LED Brightness Adjustable, 3xAAA Batteries/AC Power



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A11.1020

Student Biological Microscope






A11.1020-B Binocular Head



A11.1020-M Monocular Head



Plain Stage with Slide Clip



N.A. 0.65 with Disc Diaphragm



Upper Flexible LED, Bottom LED Brightness Adjustable



A11.1020-M



A11.1020-B

	A11.1020-M	A11.1020-B
Head	Monocular Head, 45° Inclined	Binocular Head , 45° Inclined
Eyepiece	WF10x/16mm	WF25x/9mm
		Triple Nosepiece
Nosepiece		Achromatic Objective 4x/0.10
Objective		Achromatic Objective 10x/0.25
		Achromatic Objective 40x/0.60, Spring
		Plain Stage with Slide Clip, Size 95*95mm
Stage		Coaxial Coarse & Fine Adjustment
Focusing		N.A. 0.65 with Disc Diaphragm
Condenser		Upper Flexible LED, Bottom LED Brightness Adjustable
Illumination		Power by AC 110-240V or DC 3xAAA Batteris (Batteries Not Included)
Power		Total Height 30cm, Base 19x12cm
Size		
Optional Accessories		
Vernier Scale	Sliding Ruler	A54.1002-B



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China