



## OPTO EDU A12.1005 Achromatic LED Trinocular Biological Laboratory Microscope

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1005
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



### Product Specification

- Head: Siedentopf Head
- Eyepiece: Wide Field Eyepiece WF10x/ 18
- Objective: Achromatic Objective
- Nosepiece: Backward Quadruple Nosepiece
- Working Stage: Double Layers Mechanical Stage 132x142mm/ 75x40mm
- Focusing: Coaxial Coarse & Fine Knob, Moving Range 24mm
- Condenser: Abbe Condenser N.A.1.25 With Iris Diaphragm & Filter
- Illumination: LED 3W, Brightness Adjustable
- Highlight: **OPTO-EDU Biological Laboratory Microscope, Achromatic LED Biological Laboratory Microscope**  
**Achromatic LED Trinocular Biological Laboratory Microscope**

OPTO-EDU



### More Images

OPTO-EDU



OPTO-EDU





## A12.1005 Laboratory Biological Microscope



### A12.1005 Details



Interpupillary Distance 180° Adjustable



Achromatic Objective



Double Layers Mechanical Stage





Abbe Condenser N.A.1.25



LED 3W, Brightness Adjustable

**A12.1005 Specification**



Unit: mm

A12.1005 Biological Microscope		A	B
Head	Siedentopf binocular head, inclined at 30°, interpupillary distance 48-75mm	•	
	Siedentopf trinocular head, inclined at 30°, interpupillary Distance 55-75mm		•
Eyeiece	Wide field eyepiece WF10x/ 18	•	•
Objective	Achromatic objective 4x	•	•
	Achromatic objective 10x	•	•
	Achromatic objective 40x	•	•
	Achromatic objective 100x	•	•
Nosepiece	Backward quadruple nosepiece	•	•
Stage	Double layers mechanical stage 132x142mm/ 75x40mm	•	•
Focusing	Coaxial coarse & fine knob, moving range 24mm, fine division 0.004mm	•	•
Condenser	Abbe condenser N.A.1.25 with Iris diaphragm & filter	•	•
Illumination	LED 3W, brightness adjustable	•	•
<b>Optional Accessories</b>		<b>Item No.</b>	
Eyeiece P16x/11		A51.1002-16	

Phase Contrast Kit	A5C.1009
Polarization Attachment	A5P.1003
Illumination Halogen lamp 6V/ 20W, brightness adjustable	A56.1003
Dark Field Condensor	A5D.1005



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China