



OPTO-EDU A23.3665 0.65-6.5x 1:10 Zoom Stereo Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3665
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Optical System: Greenough
- Head: Binocular / Trinocular Head, Interpupillary Distance 55~75mm
- Magnification: 0.65x~6.5x Continuous Zoom Body, Total Magnification 1.95x~325x
- Zoom Ratio: 1:10
- Eyepiece: WF10x/22mm, High Eyepoint, Diopter Adjustable, Dia.30mm
- Max Field Of View: Φ 36 Mm

OPTO-EDU



More Images

OPTO-EDU



OPTO-EDU



Product Description

Continuous Zoom 0.65x-6.5x Binocular & Trinocular Head, Zoom Ratio 1:10
 Standard 6.5x-65x, Up To 1.95x-325x With Optional Eyepieces & Auxiliary Objectives
 High Eyepoint, Dia.30mm WF10x/22mm, Optional WF15x, WF20x, WF25x
 Working Distance 110mm, Up To 30~287mm With Auxiliary Objectives
 35 Kinds Stand Optional, Details See A54.3630



OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
 Tel:+8610 88696020 Fax:+8610 88696085

A23.3665
Zoom Stereo Microscope, 0.65~6.5x, 1:10



A23.3665 Specification 



A23.3665 Zoom Stereo Microscope, 0.65-6.5x, 1:10	
Optical System	Greenough
Head	Binocular / Trinocular Head, Interpupillary Distance 55~75mm

Magnification	0.65x~6.5x Continuous Zoom Body, Total Magnification 1.95x~325x (With Optional Eyepiece & Lens)
Zoom Ratio	1:10
Eyepiece	WF10x/22mm, High Eyepoint, Diopter Adjustable, Dia.30mm
Max Field of View	Φ36 mm
Working Distance	110mm, Optional 30~287mm
Bracket Holder Hole	Inner Dia. Φ76mm
Model	A23.3665-B6T Trinocular Head, Light Split Fixed 4/6 B6 Pole Stand, Size 210*274*38mm, Pole Dia.32mm, Height 240mm, No Light Source, 0.5x C-Mount
	A23.3665-B8LT Trinocular Head, Light Split Fixed 4/6, B8L Pole Stand, Size 210*274*38mm, Pole Dia.32mm, Height 240mm, Up/Bottom 3W LED, 0.5x C-Mount
	A23.3665-B8LST Trinocular Head, Light Split Fixed 4/6, B8LS Pole Stand, Size 210*274*38mm, Pole Dia.32mm, Height 240mm, Up Dual Pipe LED Brightness Separately Adjustable, Bottom 3W LED, 0.5x C-Mount
	A23.3665-B3T Trinocular Head, Light Split Fixed 4/6, B3 Pole Stand, Size 380*260*20mm, Pole Dia.32mm, Height 240mm, No Light Source, 0.5x C-Mount
Optional	Eyepiece WF15x, WF20x, WF25x, A51.3622
	Auxiliary Lens 0.5x, 0.7x, 1.5x, 2.0x, A52.3622
	C-Mount 0.35x, 1.0x, A55.3623
	35 Kinds of Stands & Brackets, See A54.3630

A23.3665 Optical Data Sheet												
Eyepiece	Auxiliary Objective Lens											
	0.3x		0.5x		0.7x		1.0x (Standard)		1.5x		2.0x	
	Working Distance 287mm		Working Distance 165mm		Working Distance 120mm		Working Distance 110mm		Working Distance 45mm		Working Distance 30mm	
	Mag.	View(mm)	Mag.	View(mm)	Mag.	View(mm)	Mag.	View(mm)	Mag.	View(mm)	Mag.	View(mm)
WF10x/22mm	1.95x	107.3	3.25x	64.4	4.55x	46	6.5x	32.2	9.75x	21.5	13x	16.1
	19.5x	11.3	32.5x	6.8	45.5x	11.2	65x	3.4	97.5x	2.3	130x	1.7
WF15x/16mm	2.92x	78.03	4.87x	46.83	6.82x	33.45	9.75x	23.42	14.62x	15.64	19.5x	11.71
	29.2x	8.21	48.75x	4.95	68.25x	8.15	97.5x	2.47	146.25x	1.67	195x	1.24
WF20x/12mm	3.9x	58.52	6.5x	35.13	9.1x	25.09	13x	17.56	19.5x	11.73	26x	8.78
	39x	6.16	65x	3.71	91x	6.11	130x	1.85	195x	1.25	260x	0.93
WF25x/9mm	4.87x	43.89	8.12x	26.34	11.37x	18.82	16.25x	13.17	24.37x	8.79	32.5x	6.59
	48.75x	4.62	81.25x	2.78	113.75x	4.58	162.5x	1.39	243.75x	0.94	325x	0.69





A23.3665-B3T



A23.3665-B3T Work With A59.5104 WIFI Digital Eyepiece Camera



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China