



Monocular Cordless Biological Microscope Mechanical Working Stage A11.0020

Our Product Introduction

Basic Information

- Place of Origin: Korea
- Brand Name: CNOEC OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.0020
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Student
- Head: Monocular
- Nosepiece: Triple
- Illumination: LED Lamp
- Power: Cordless
- Condensor: N.A. 0.65
- Highlight: **binocular biological microscope,
biological compound microscope**



for more products please visit us on cnoec.com

Product Description

Monocular Cordless Biological Microscope Mechanical Working Stage A11.0020

A11.0020 Cordless Student Microscope, newly developed in 2015, with modern design and high optical quality, is ideal tools in educational field.

- Monocular Head Rotatable
- Triple Nosepiece
- Clip Working Stage
- N.A 1.25 Condenser, Disc Diaphragm
- Colorful Focusing Knob, Light Collector and Power Switch
- LED Power On/Off Indicator
- LED Brightness Adjustment Indicator
- Cordless Rechargeable Battery Built-in

A11.0020 Cordless Student Microscope		A11.0020	-60	-100
Head	Monocular Head, 30° Inclined, 360° Rotatable			
Eyeiece	WF10x/18mm Wide Field Eyepiece			
Nosepiece	Triple Nosepiece, Ball Bearing Inner Locating			
Objective	Achromatic Objective 4x, 10x, 40x (S)	o		
	Achromatic Objective 10x, 40x (S), 60x (S)		o	
	Achromatic Objective 10x, 40x (S), 100x (S)			o
Working Stage	Single Layer Stage 125*120mm, With Specimen Clips			
Focusing	Coaxial Coarse & Fine Focus Knob, With Tension Adjustable And Limit Stopper, Minimum Division 0.002mm			
Condenser	N.A. 0.65 Condenser, 6 Holes Disc Diaphragm			
Light Source	High Light LED, Brightness Adjustable, With LED Indicator			
Power	Cordless Built-in Rechargeable Battery 100~240V Wide Range Voltage Power On/Off LED Indicator			
Packing	Net Weight 4.28 Kgs/pc, Packing 2 pcs/Carton, G.W. 10Kgs, Size 34*48*43cm			



High Quality Ach. Objective 4x,10x,40x,60x,100x Available



Colorful Focusing Knob & Power Switch, Power Indicator




Modern Designed Brightness LED Indicator



Round Edge Working Stage Protect Hand From Hurt

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China