



Polarized Light Microscope Trinocular Head Round Stage A15.0013

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE
- Model Number: A15.0013
- Minimum Order Quantity: 1
- Price: Negotiation
- Packaging Details: Carton
- Delivery Time: 1 month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000



Product Specification

- Application: Educational
- Certification: CE
- Head: Trinocular
- Illumination: LED Lamp
- Nosepiece: Triple
- Stage: Round Stage
- Highlight: **student compound microscope ,
highest magnification microscope**

for more products please visit us on cnoec.com

Product Description

Polarized Light Microscope Trinocular Head Round Stage A15.0013

A15.0013 Polarizing Microscope, with modern design and high optical quality, is ideal tools in educational field.

- Trinocular Head
- Triple Nosepiece
- Round Double Layer Mechanical Working Stage
- N.A.1.25 Condenser, Iris Diaphragm
- Polarizer & Analyzer, Built-in Bertrand Lens, Compensator
- Colorful Focusing Knob, Light Collector and Power Switch
- LED Power On/Off Indicator
- LED Brightness Adjustment Indicator
- Cordless Rechargeable Battery Built-in

Specification

A15.0013 Polarizing Microscope	
Head	Trinocular Head, 30° Inclined, 360° Rotatable
Eyepiece	WF10x/18mm, WF15x/13mm
Nosepiece	Triple Nosepiece, Ball Bearing Inner Locating
Objective	Achromatic Objective 4x, 10x, 40x (S)
Working Stage	Round, Diameter 142mm, Rotatable 360°, Division 1°, Vernier division 6, Lockable
Focusing	Coaxial Coarse & Fine Focus Knob, With Tension Adjustable And Limit Stopper, Minimum Division 0.002mm
Condenser	N.A. 1.25 Condenser, Iris Diaphragm, Filter
Polarizing	Analyzer : 0-90° Rotatable, Can Be Moved Out Of Optical Path For Single Polarizing Observe Bertrand Lens Built-in Compensator : λ Slip (First Class Red) , $1/4\lambda$ Slip, Quartz Wedge Plate Polarizer : Mounted Under Abbe Condenser, 360° Rotatable, With Scale, Lockable
Light Source	High Light LED, Brightness Adjustable, With LED Indicator
Power	Cordless Built-in Rechargeable Battery 100~240V Wide Range Voltage Power On/Off LED Indicator
Packing	Net Weight 4.80Kgs/pc, Packing 2 pcs/Carton, G.W. 10Kgs, Size 34*48*43cm



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627
 sale@optoedu.com
cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China