



LED Light Source Compound Microscope Parts Double Pipe Light , Dual Knob A56.2412

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.2412
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Light Source: LED
- Power: 3W LED
- Brightness: Adjustable
- Highlight: **microscope accessory,
compound microscope parts**

for more products please visit us on cnoec.com

Product Description

LED Light Source Compound Microscope Parts Double Pipe Light , Dual Knob A56.2412


- Double 3W LED Pipe
- Brightness Adjustable

A56.5202-4 Four Color LED Ring Light	
Base	6*11*14cm
Pipe Length	55cm
Pipe Connection	Fixed or Removeable
Power	3W
Input Voltage	Wide Range Voltage100-240V
Packing	18*18*30cm, 2.5Kgs



A56.2412





Input 5V/2A or Input 100V-240V is available



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China