



## Trinocular 5.0M Digital Optical Microscope USB LED Light A31.1535

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A31.1535
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Digital Microscope
- Certification: CE
- Resolution: 5.0M
- Illumination: LED
- Output: USB2.0
- Sensor: COMS
- Highlight: **digital usb microscope,**  
**portable digital microscope**



### More Images



## Product Description

Trinocular 5.0M Digital Optical Microscope USB LED Light A31.1535

- USB2.0 Output
- 5.0M resolution
- Led Light Source


### A31.1535 Digital Biological Microscope

Camera	1/2.5"CMOS , Active 5.0M, Resolution 2592 x 1944, USB2.0 port
Head	Trinocular, Inclined at 30°, Interpupillary 48-75mm
Eyeiece	WF10x/18mm
Nosepiece	Backward Quadruple Nosepiece
Objective	Achromatic Objectives 4x, 10x, 40x, 100x
Condensor	Abbe Condensor NA1.25 with Diaphragm , Filter
Focusing	Coaxial Coarse and Fine Adjustment, Moving Range 22mm
Working Stage	Double Layer Mechanical Stage 140x130mm
Light Source	LED Light, Brightness Adjutable
Adapter	0.5x C-mount



A31.1535 Digital Biological Microscope, 5.0M







**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China