



## Student Biological Compound Microscope LED Light Source A11.6112

Our Product Introduction

### Basic Information

- Place of Origin: Korea
- Brand Name: OSUN
- Certification: CE, Rohs
- Model Number: A11.6112
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Application: Student
- Head: Monocular
- Nosepiece: Triple
- Illumination: LED Lamp
- Focus System: Coaxial Coarse And Fine Focus
- Condensor: N.A. 1.25
- Highlight: **high powered microscope ,  
binocular biological microscope**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Student Biological Compound Microscope LED Light Source A11.6112

**A11.6112 Student Microscope**, newly developed in 2016, with modern design and high optical quality, is ideal tools in educational field.

- Quaduple Nosepiece
- Double Layer Mechanical Stage
- N.A 1.25 Condenser, Disc Diaphragm
- LED Brightness Adjustment Indicator

A11.6112 Biological Microscope		B	T
Head	Binocular	o	
	Trinocular		o
Eyeiece	WF10x/18		
Nosepiece	Backward Quaduple Nosepiece		
Objective	Achromatic Objectives 4x, 10x, 40x & 100x (oil)		
Stage	Double Layer Mechanical Stage		
Focus System	Coaxial Coarse And Fine Focus		
Condenser	Abbe Condenser 1.25, Center Adjustable, Height Adjustable		
Illumination	LED 3W Brightness Adjustable, Wide Range Power Adapter		



A11.6112



Side View



Front View






**Opto-Edu (Beijing) Co., Ltd.**

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China