



Abbe Refractometer Microscope Accessories Binocular B31.0101

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: B31.0101
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Binocular
- Model: Protoble
- Display: Digital
- Highlight: compound microscope parts, parts of the microscope



for more products please visit us on cnoec.com

Product Description

Abbe Refractometer Microscope Accessories Binocular B31.0101

- Monocular and Binocular
- Portable

| Specifications | B31.0101-M | B31.0101-B |
|---|----------------|----------------|
| Head | Monocular Head | Binocular Head |
| Refractive index ND measurement range | 1.300 1.700 | |
| Refractive index ND measurement accuracy | 0.0002 | 0.0003 |
| Refractive index ND Min. division | 0.001 | 0.0005 |
| Glucose concentration (%) measurement range | 0 95 | |
| Glucose concentration (%) Min. division | 0.5 | 0.25 |



B31.0101-B
Abbe Refractometer

www.cnoec.com



B31.0101-M
Abbe Refractometer



www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627

sale@optoedu.com

cnoec.com

