



Infinity Plan Compound Optical Microscope with Halogen Lamp A12.2701-AP

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.2701- AP
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Type: Biological Microscope
- Optical System: Infinity Correct
- Light Source: Halogen Lamp
- Magnification: 1000x
- Application: School
- Head: Binocular
- Highlight: optical compound microscope, wide field microscope



More Images



for more products please visit us on cnoec.com

Product Description

Infinity Plan Compound Optical Microscope with Halogen Lamp A12.2701- AP

A12.2701- AP Laboratory Biological Microscope, Brand New Developed in 2016, is the best choice for your lab and school!

Perfect Optical System

Perfect Appearance Design

Perfect Ergonomics Design

Infinity Correct Optical System Supply High Resolution And Perfect Definition

Specification	A12.2701-AP
Head	Binocular Head 30° Inclined, 360° Rotatable, Interpupillary Distance 48-75mm, Diopter Adjustment -/+5D,
Eyepiece	WF10X/18mm, High Eyepoint, Eye Cup,
Objective	Infinity Plan Achromatic Objectives 4x/0.13, 10x/0.30, 40x/0.70 (Spring), 100x/1.25 (Spring, Oil)
Focus System	Low Positioning, Coaxial Coarse/Fine Focus System, Coarse Knob Tension Adjustable, 14mm Per Rotation Fine Knob 0.1mm Per Rotation, Minimum Reading 1 Microns On Left-Side Fine Knob
Nosepiece	Quadruple Inward Nosepiece
Stage	Double Layer Mechanical, Size 156x138mm, Moving Range 76x54mm, Low Positioned X-Y Coaxial Control Knob, Can Hold 2 Specimen Slides
Condenser	Abbe N.A. 1.25 With Aperture Iris Diaphragm And With Socket For phase Contrast And Dark Field Slider
Filter	Blue Filter
Light Source	Built-In Transmitted, Kohler, 6V 20W Halogen, With Brightness Adjustable, Independent Reflected Illuminator, With Brightness Adjustable And Rechargeable
Spare Parts	Wooden Case With Lock And Keys, Immersion Oil 100ml, Spare Bulb 2pcs, Cleaning Paper 10*15cm(50pcs/Book) 2 Books, Cleaning Brush 2pcs
Power Supply	EU Plug, 220V 50Hz




High precision Quadruple Nosepiece (Quintuple are optional)



Double Layer Mechanical, Size 156x138mm, Moving Range 76x54mm. Low Position, Right Hand Control, X-Y Coaxial Control Knob, Two Slides Holder



Abbe Condenser N.A. 1.25



Bright Field



Dark Field



Phase Contrast



Polarization



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China