



Student Biological Microscope With S-LED Illumination A11.1001

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1001
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Product Name: Microscope
- Theory: Biological Microscope
- Objective: 185 Achromatic Objective
- Applications: Student; Laboratory
- Nosepiece: Triple & Quadruple Nosepiece
- Magnification: 40x-1000x
- Highlight: **high powered microscope ,
binocular biological microscope**



More Images



Product Description

Specifications

| A11.1001 Biological Microscope | | A | B | C |
|--------------------------------|--|-----------------|---|---|
| Head | Monocular Head, Inclined 45° | ● | | ● |
| | Binocular Head, Inclined 30° | | ● | |
| Eyepiece | Wide Field Eyepiece WF10x/18mm | ● | ● | ● |
| Objective | 185 Achromatic Objective 4x | ● | ● | ● |
| | 185 Achromatic Objective 10x | ● | ● | ● |
| | 185 Achromatic Objective 40x | ● | ● | ● |
| | 185 Achromatic Objective 100x | | ● | ● |
| Nosepiece | Triple Nosepiece | ● | | |
| | Quadruple Nosepiece | | ● | ● |
| Focusing | Coarse Focusing | ● | | |
| | Separate Coarse & Fine Focusing | | ● | ● |
| Stage | Plain Stage With Slide Clips, Stage Size 125*120mm | ● | | |
| | Attachable Mechanical Stage Size 110*120mm/ 60*30mm | | ● | ● |
| Condenser | Single Lens With Iris | ● | | |
| | Abbe N.A.1.2 With Iris, Spiral Up & Down | | ● | ● |
| Illumination | S-LED Illumination, Brightness Adjustable | ● | ● | ● |
| Filter | Blue Filter | | ● | ● |
| | Green Filter | | ● | ● |
| Optional Accessories | | Item No. | | |
| Head | Dual-View Head | A53.1001-1 | | |
| Eyepiece | Wide Field Eyepiece P16*11 | A51.1002-16 | | |
| Illumination | Plan-concave Mirror | A56.1001-1 | | |

A11.1001 Biological Microscope



A11.1001-A

A11.1001-C

A11.1001-B



Opto-Edu (Beijing) Co., Ltd.

