



110x115x50mm Microscope Accessories Magnet Stand A54.2501 For A13.2501

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A54.2501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope Accessories
- Certification: CE|Rohs
- Product Name: Microscope Magnet Stand
- Work Table Size: 90x90mm
- Inner Diam: 65/50mm
- Size: 110x115x50mm
- Highlight: **compound microscope parts, parts of the microscope**



Product Description

110x115x50mm Microscope Accessories Magnet Stand A54.2501 For A13.2501

This stand is suitable to fix onto A13.2501 microscope, it is flexible to move and operate.

Specification

A54.2501 Magnet Stand For A13.2501	
Work Table size	90x90mm
Inner Diam	65/50mm
Size	110x115x50mm
Used	Attach To Metal Surface



A54.2501 Magnet Stand



A54.2501 Magnet Stand



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China