



A56.2618 Microscope Accessories Cool Light Fiber Optic Halogen Illuminator

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.2618
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Model: A56.2618
- Power: 21V 150W With Philips Bulb 13629
- Product Name: Cool Light Fiber Optic Halogen Illuminator
- Packing Size: 285x230x255mm
- Size: 230x101x150mm
- Filter Life: Light Blue (Optional: Navy Blue, Red, Green)
- Highlight: **microscope accessory, parts of the microscope**

for more products please visit us on cnoec.com

Product Description

A56.2618 Microscope Accessories Cool Light Fiber Optic Halogen Illuminator


- For observing temperature-sensitive specimens - Gooseneck lights can be positioned independently

Specification

A56.2618 Cool Light Fiber Optic Halogen Illuminator		
Power	21V 150W with Philips Bulb 13629	
Input voltage	76V-265V, 50Hz (110V optional)	
CT	3000K	
Illuminance	Approximately 100000Lx	
Filter	Light Blue (Optional: Navy Blue, Red, Green)	
Fiber Optic Interfacitive	16mm (5/8")	
Size	230x101x150mm	
Net Weight	2.1kgs	
Gross Weight	2.7kgs	
Packing Size	285x230x255mm	
Optional Accessories		Item No.
Fiber Optic Light Guide (Gooseneck type)	Metal-sheathed Fiber Optic Bifurcated Light Guides With Condenser, Fiber Length 550mm, Fiber Diameter 8mm, Standard Interface 5/8"	A56.2610-03
	Metal-sheathed Fiber Optic Single Light Guide with Condenser, Fiber Length 550mm, Fiber Diameter 8mm, Standard Interface 5/8"	A56.2610-04
	Metal-sheathed Fiber Optic Ring Light Guide With Condenser, Fiber Length 550mm, Fiber Diameter 8mm, Standard Interface 5/8"	A56.2610-02



A56.2618
Cool Light Fiber Optic
Halogen Illuminator





 0086 13911110627  sale@optoedu.com  cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China