



Finite Optical System Biological Microscope With Coaxial Coarse / Fine Focusing , A12.1035

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: A12.1035
- Minimum Order Quantity: 1
- Price: Negotiation
- Packaging Details: Carton
- Delivery Time: 1month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000



Product Specification

- Product Name: Biological Microscope
- Optical System: Finite Optical System
- Stage: Double Layers Mechanical Stage
- Focusing: Coaxial Coarse And Fine Adjustment
- Condenser: Abbe Condenser N.A.1.25
- Illumination: S-LED Illumination
- Highlight: **high powered microscope , biological compound microscope**



for more products please visit us on cnoec.com

Product Description

Finite Optical System Biological Microscope With Coaxial Coarse and Fine Focusing A12.1035

Specification

A12.1035 Biological Microscope		B	T	B E	TE
Optical System	Finite Optical System				
Eyepiece	WF10X/18				
Viewing Head	Monocular Head Inclined at 30°	○	○	○	○
	Seidentopf Binocular Viewing Head Inclined at 30°, Interpupillary 47-78mm	●		●	
	Seidentopf Trinocular Viewing Head Inclined at 30°, Interpupillary 47-78mm		●		●
	Digital Viewing Head Inclined at 30°, Interpupillary 47-78mm, 5M pixels, Support Wifi	○	○	○	○
	Digital Viewing Head Inclined at 30°, Interpupillary 47-78mm, 5M pixels, 8" LCD. Android OS	○	○	○	○
Achromatic Objective	4X N.A=0.1 WD=37.5	●	●	●	●
	10X N.A=0.25 WD=7.316	●	●	●	●
	20X N.A=0.4 WD=13	○	○	○	○
	40X N.A=0.65 WD=0.632	●	●	●	●
	100X N.A=1.25 WD=0.198	●	●	●	●
Nosepiece	Backward Quadruple Nosepiece	●	●	●	●
Stage	Double Layers Mechanical Stage 140X140mm, Moving Range 75X50mm				
Condenser	NA1.25 Abbe Condenser	●	●	●	●
	Dark Field Condenser (Dry, Oil)	○	○	○	○
Focusing System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Moving Range 25mm				
Illumination	Halogen Lamp 6V/20W	○	○	○	○
	S-LED Illumination	●	●	●	●
	S-LED Illumination with Recharger	○	○	○	○
	Kohler Illumination	○	○	○	○
	ECO Function (Auto Power Off After 30 Mins No Operation)	○	○	●	●
	Plane-concave Mirror	○	○	○	○
Fluorescent Attachment	Fluorescent Illumination, LED Lamp	○	○	○	○
	Fluorescent Illumination, 100W HBOLamp	○	○	○	○
Polarzation Set		○	○	○	○
Phase Contrast Kit		○	○	○	○

Note: ● Standar Outfit , ○ Optional

www.cnoec.com



ECO Function Newly Developed!
Automatic Power Off When
The People Keep Away Over 30Mins

A12.1035-BE Microscope With ECO Function

www.cnoec.com



A12.1035 Microscope Back View



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

