



## 400X 360° Rotatable Trinocular Inverted Optical Microscope Biological Microscope A14.0912

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A14.0912
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Trinocular
- Illumination: Halogen Lamp
- Magnification: 400X
- Product Name: Inverted Biological Microscope
- Highlight: **inverted biological microscope, inverted light microscope**

## Product Description

### 400X 360°Rotatable Trinocular Inverted Biological Microscope A14.0912

#### Specification

A14.0912 Inverted Microscope		
Optical System	Infinity Color Corrected Optical System	
Head	Trinocular Head, Split Ratio 50:50, 45° Inclined, 360° Rotatable, IPD 50-75mm,	
Eyepiece	Plan WF10x/22mm, Diopter Adjustable, High Eye-point,	
Objective	Infinity LWD Plan Achromatic Objective 4x/10x/40x Infinity LWD Positive Phase Contrast Objective 20x	
Focus System	Coaxial Focus System With Tension Adjustment, Coarse Range: 9mm, Fine Precision 2μm,	
Nosepiece	Quadruple Nosepiece	
Condenser	LWD N.A. 0.30, W.D. 72mm, Dismountable, Phase Contrast Slide, Centering Telescope	
Stage	Right Hand Controls, Low Positioned, Fixed Mechanical stage, Size 215x250mm, φ35 Petri Dish Holder, Metal Plate And Glass Plate	
Light Source	LED Lamp, 5W	
	Light Indicator	
	Brightness Adjustable	
	ECO Infrared Induction System (Auto on/off)	
Filter	External Wide Voltage Transformer, 100V-240V Input, 12V5A Output,	
	LBD color temperature change filter(φ45mm) IF550 green interference filter(φ45mm)	
Optional Accessories		Item No.
Eyepiece	Plan WF15x/16mm, Diopter Adjustable, High Eye-point	
Stage	Mechanical mobile ruler, moving range120(X)*78(Y)mm Stage extension plate	A54.0945-03
Adapter	035x C-Mount	
	0.5x C-Mount	
	0.65x C-Mount	
	1x C-Mount	



A14.0912 Inverted Biological Microscope



**ECO Infrared Induction System,Auto on/off For Energy Saving**



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China