



Monocular WF10x Student Biological Microscope with N.A.0.65 Condenser A11.6111

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.6111
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student|School
- Certification: CE|Rohs
- Drawtube: Monocular
- Illumination: LED
- Magnification: 100x-400x
- Product Name: Biological Microscope
- Highlight: **high powered microscope ,
biological compound microscope**



Product Description

Monocular WF10x Student Biological Microscope with N.A.0.65 Condenser A11.6111

Specification

A11.6111: Monocular, WF10x, Triple, 4x10x40x, Clip Stage, N.A.0.65 Condenser, 6 Holes Disc Diaphragm, Coaxial Coarse & Fine Focusing, LED, With Packing Bucket

A11.6111-2: Monocular, WF10x, WF25x, Triple, 4x10x40x, Clip Stage 90x90mm, N.A.0.65 Condenser, 6 Holes Disc Diaphragm, Coaxial Coarse & Fine Focusing, LED or 3pcs AA batteries

A11.6111-3: Monocular, WF10x, WF25x, Triple, 4x10x40x, Clip Stage 90x90mm, N.A.0.65 Condenser, 6 Holes Disc Diaphragm, Coaxial Coarse & Fine Focusing, (Blue Plain Knob), LED or 3pcs AA batteries

A11.6111-4: Monocular, WF10x, WF25x, Triple, 4x10x40x, Clip Stage 90x90mm, N.A.0.65 Condenser, 6 Holes Disc Diaphragm, Coaxial Coarse & Fine Focusing, (Blue Plain Knob), LED or 3pcs AA batteries (rechargeable)



A11.6111 Student Microscope



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China