



## Finity Optical System WF10X/18mm Binocular Microscope A11.1550 With Halogm Lamp Illumination

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1550
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Optical System: Finity Optical System
- Certification: CE|Rohs
- Focusing System: Coaxial Coarse And Fine Adjustment
- Illumination: Halogm Lamp
- Magnification: 1000X
- Product Name: Biological Microscope
- Highlight: **binocular biological microscope, biological compound microscope**



### More Images



## Product Description

### Finity Optical System WF10X/18mm Biological Microscope A11.1550 With Halogm Lamp Illumination

#### Specifications:

A11.1550 Biological Microscope		
	A11.1550-B	A11.1550-T
Optical System	Finity Optical System	
Head	Seidentopf Binocular Viewing Head, Inclined at 30° ,360°Rotatable	Seidentopf Trinocular Head,Inclined at 30° ,360°Rotatable
Eyeiece	Wide Field Eyeiece WF10X/18mm	
Objective	Achromatic Objectives 4X 10X 40Xs 100Xs(oil)	
Noseiece	Quadruple Noseiece	
Condenser	Abbe Condenser N.A.1.25 with Diaphragm, Filter	
Focusing System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke36mm	
Stage	Double Layer Mechanical Stage 142x132mm , Moving Range75mmx50mm	
Illumination	3V/20W Halogm Lamp,Brightness Adjustment	
Adapter	N/A	Eyeiece Tube Adapter
A11.1550 Biological Microscope Optional Accessories		
Eyeiece	Wide Field Eyeiece WF16X/13mm	A51.1503
	Wide Field Eyeiece WF20X/11mm	A51.1503
Objective	Plan Achromatic CF Objectives 4X 10X 40Xs 100Xs(oil)	A52.1510
Noseiece	Quintuple Noseiece	A54.1502-5A
Illumination	3W/LED Adjustment	A56.1503
	Rechargeable Battery	A55.1550-10
CCD Adapter	1.0x	A55.1550-10
	0.5x	A55.1550-05



A11.1550-B Biological Microscope



**A11.1550-T Biological Microscope**



**A11.1550-T + Digital Camera**



**A33.1550-A**  
A11.1550-T + LCD Digital Camera



**A33.1550-B**  
A11.1550-T + LCD Digital Camera



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China