



## A11.1323 WF10x Biological Microscope With Single Lens N.A.0.65 Condenser

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1323
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Student|School
- Stage: Plain Stage With Specimen Clips
- Drawtube: Monocular
- Illumination: LED
- Magnification: 100x-400x
- Product Name: Biological Microscope
- Highlight: **high powered microscope ,  
binocular biological microscope**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Monocular WF10x Student Biological Microscope With Single Lens N.A.0.65 Condenser A11.1323

#### Introduction


A11.1323 is a very good choice for student, especially appreciated when budget is limited. All models are equipped with AC110/220V and internal battery, which can meet your different needs. Achromatic 4X, 10X, 40X objective enable utmost 400X optical magnification, which can clearly observe biological slices, biological cells, bacteria and living tissue culture etc.



#### Parameter

A11.1323 Student Microscope	
Head	Monocular Head, 45° Inclined, 360° Rotatable
Eyepiece	WF10X
Objective	Achromatic 4X
	Achromatic 10X
	Achromatic 40X, Spring
Nosepiece	Triple Revolving Nosepiece
Stage	Plain Stage With Specimen Clips, Size 90x90mm
Condenser	Single Lens N.A.0.65 With Six-holes Disc Diaphragm
Focusing	Coaxial Coarse & Fine Focusing Adjustment, With Safty Stop
Light Source	LED Light
Power Supply	External AC/DC Adaptor, Or By 3xAA Batteries(Not Included)
Size	16x11x31cm, Net Weight 1.7KG

#### Product Display



### Product Details




## 01 MONOCULAR HEAD

Rotatable 360° with WF10X eyepiece.

## 02 NOSEPIECE

Quadruple nosepiece with precise click stop, equipped with achromatic objective 4X, 10X, 40Xs.



## 03 WORKING STAGE

Plain stage with clip and single lens.

## 04 DISC DIAPHRAGM

With single lens condenser, easy to adjust the light source.



## 05 FOCUSING

Coaxial coarse & fine focusing adjustment.


## 06 POWER SUPPLY

It is powered by outlet 110-240V adapter, or by 3 pcs AA batteries for cordless usage.



### Product Application

Student microscope is designed for tertiary institutions and primary and secondary schools, which convenient to observe biological slices, biological cells, bacteria and living tissue culture, fluid precipitation observation and research. It is widely used in teaching demonstration, biochemical experiments, clinical research ect.



### Packaging & Shipping

	1 Set/Carton	10 Sets/Carton
Package	34x23x16CM	58x49x36CM
Gross Weight	1.7KG	17KG
Packing	Polyfoam Protection And Strong Carton.	
Shipping Time	5~20 Days After Receive The Deposit.	



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China