



Opto Edu Ent Operating Microscope A41.1901-A Eyetube Diopter Adjustable 6x

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A41.1901-A
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: 17.24x
- Head Magnification: 6x
- Interpupillary Distance: 50~80mm Adjustable
- Diopter: Both Left & Right Eyetube Diopter Adjustable $\pm 5D$
- Zoom Body: 3 Steps Zoom: 0.6x, 1.0x, 1.6x
- Focusing: Fine Focusing Range 10mm
- Highlight: **zoom ent operating microscope,**
6x ent operating microscope

for more products please visit us on cnoec.com

Product Description

OPTO-EDU A41.1901-A Eyetube Diopter Adjustable Operating Microscope For ENT Surgery 6x

A41.1901 Series Operation Microscope Specification

Straight Binocular 6x, Zoom Body 0.6x, 1.0x, 1.6x, Objective F=200mm, 300mm, Total Magnification 3x~12x, Built-in Blue, Orange Filter, Fine Focusing Range 10mm, 10W LED, Column + Five Star Base With Caster

 OPTO-EDU (BEIJING) CO., LTD. 

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A41.1901

Series Operating Microscope,
Dental, ENT, Ophthalmology, Gynecology





Model	A41.1901-A	A41.1901-C	A41.1901-CZ	A41.1901-B	A41.1901-D
Binocular Head	Straight			45° Tilted	
Main Usage	ENT, Surgery			Dental, Ophthalmology, Gynecology	
Total Magnification	With Objective F=200mm: 4.7x, 7.5x, 12x With Objective F=300mm: 3.0x, 5.0x, 8.0x				
Total Magnification (Optional)	With Objective F=250mm: 3.7x, 6.0x, 9.6x With Objective F=400mm: 2.3x, 3.8x, 6.0x				
Eyepiece	17.24x				
Head Magnification	6x				
Interpupillary Distance	50~80mm Adjustable				
Diopter	Both Left & Right Eyetube Diopter Adjustable $\pm 5D$				
Zoom Body	3 Steps Zoom: 0.6x, 1.0x, 1.6x				
Objective	F=200mm For Ear, Nose F=300mm For Ear, Nose			F=200mm For Ophthalmology F=300mm For Gynecology	
Objective (Optional)	F=400mm For Throat			F=250mm For Dental	
Focusing	Fine Focusing Range 10mm				
Filter	Green, Yellow Filters Built-in				
Balancing Arm	2-Part Arm With Universal Joints				
Stand	2-Part Column + Five Star Base With Caster	Table Mount Clamp	Column Mount, Customized To Customer	2-Part Column + Five Star Base With Caster	Table Mount Clamp
Illumination System	Coaxial Illumination, 10W LED, Brightness Adjustable, >30000 Lux				
Power Supply	Wide Range AC100~240V, DC12V 3A, 36W				
Video Adapter	Not Available				
A41.1901 Series Operation Microscope Optional Accessories					
Objective	F=250mm, For A41.1901-B, -D			A52.1901-250	
	F=400mm, For A41.1901-A, -C, -CZ			A52.1901-400	
Mount	Table Mount Clamp			A54.1901-T	
	Wall Mount Bracket			A54.1901-W	
Monitor Mount	Monitor Mount: 1 Fix Arm And 1 Gas Spring Arm With Tension Adjustable. Band Clamp. Fit 15-27" Screen. VESA Standard 75x75xmm & 100x100mm. Load Capacity 0~8 Kgs, Swivel 180°, Tilt -80°~+90°, Rotation 360°				A54.1910-112B
	Monitor Mount: 1 Fix Arm And 1 Mechanic Spring Arm With Tension Adjustable. Band Clamp. Fit 15-27" Screen. VESA Standard 75x75xmm & 100x100mm. Load Capacity 0~8 Kgs, Swivel 180°, Tilt -80°~+90°, Rotation 360°				A54.1910-312B

A41.1901 Series Operation Microscope Optical Characteristics												
Focal Length (mm)	F=200mm			F=250mm (Optional)			F=300mm			F=400mm (Optional)		
Zoom Body	1.6x	1.0x	0.6x	1.6x	1.0x	0.6x	1.6x	1.0x	0.6x	1.6x	1.0x	0.6x
Total Magnification	12x	7.5x	4.7x	9.6x	6.0x	3.7x	8.0x	5.0x	3.0x	6.0x	3.8x	2.3x
Linear Field (mm)	15.8	25.3	40.6	19.7	31.6	50.7	23.6	37.9	60.8	31.5	50.5	81.0
Exit Pupil Dia. (mm)	1.04	1.66	1.66	1.04	1.66	1.66	1.04	1.66	1.66	1.04	1.66	1.66
Exit Pupil Distance (mm)	15.30											


Resolution (LP/mm)	67.0	44.5	29.7	60.0	35.4	27.0	47.2	31.5	23.6	35.0	27.0	19.8
-----------------------	------	------	------	------	------	------	------	------	------	------	------	------



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China